

Withdrawn  
Provided for Historical Reference only

Hydrogen Sulfide

---

Method Number	ID-141
Matrix	Air
OSHA Permissible Exposure Limits Hydrogen Sulfide (Final Rule Limits):	10 ppm (TWA) 15 ppm (STEL)
Hydrogen Sulfide (Transitional Limits):	20 ppm (Ceiling) 50 ppm (Peak 10 min exposure)
Collection Procedure	A calibrated personal sampling pump is used to draw air through a filter impregnated with silver nitrate which converts the hydrogen sulfide to silver sulfide.
Recommended Air Volume	2 to 6 L
Recommended Sampling Rate Peak, Ceiling or STEL samples:	0.2 L /min
TWA samples	0.1 L/min
Analytical Procedure	The sample is analyzed as sulfide by differential pulse polarography.
Detection Limits Qualitative:	0.4 ppm ( 2L air volume)
Quantitative:	0.9 ppm (2L air volume)
Precision and Accuracy Validation Level:	10.1 to 39.2 ppm
CV <sub>T</sub>	0.038
Bias	-0.031
Overall error	±10.7%
Method Classification	Validated Method
Chemist	Ted Wilczek
Date (Date Revised)	1983 (December 1989)

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
Provided for Historical Reference only

---

---

---

Commercial manufacturers and products mentioned in this method are for descriptive use only and do not constitute endorsements by USDOL-OSHA. Similar products from other sources can be substituted.

---

---

---

Branch of Inorganic Methods Development  
OSHA Technical Center,  
Sandy, UT- 84070

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
Provided for Historical Reference only

## 1. Introduction

This method describes the collection of airborne hydrogen sulfide ( $H_2S$ ) in the workplace on a silver nitrate ( $AgNO_3$ )-impregnated filter and analysis by differential pulse polarography (DPP).

### 1.1. History

Previously,  $H_2S$  was collected in a midjet impinger containing an alkaline suspension of cadmium hydroxide. The sulfide was precipitated as cadmium sulfide ( $CdS$ ) and subsequently analyzed by the methylene blue calorimetric procedure (8-1.). Due to the photosensitivity of  $CdS$ , it was necessary to protect the impinger sample from light at all times. Also, the impinger base and stem could contain  $CdS$  deposits after sample collection. This required the entire impinger sample be sent to the lab for analysis. Problems associated with impinger sampling were additional incentives to develop a more acceptable sampling method.

### 1.2. Principle

Hydrogen sulfide is collected on a Whatman 4 filter paper (Whatman Labsales, Hillsboro, OR) which has been impregnated with  $AgNO_3$ . The  $H_2S$  reacts with the  $AgNO_3$  to form silver sulfide, a greyish-black precipitate (8.2., 8.3.). The silver sulfide is dissolved in an alkaline cyanide solution and analyzed for sulfide by DPP using a dropping mercury electrode (DME).

### 1.3. Advantages and Disadvantages

- 1.3.1. The sampling device is small, portable, and involves no liquids. Humidity does not significantly affect the sampling efficiency of the device.
- 1.3.2. The sulfide formed is stable and non-volatile. Desorption and preparation of samples for analysis involve simple procedures.
- 1.3.3. Collected samples are analyzed by means of a quick analytical method.
- 1.3.4. This method has adequate sensitivity for measuring workplace atmospheric concentrations of  $H_2S$  for Short-Term Exposure Limit (STEL), Ceiling, or Time Weighted Average (TWA) measurements.
- 1.3.5. The analysis is specific for sulfide in the presence of other organic or inorganic sulfur compounds.
- 1.3.6. One disadvantage of the method is the amount of  $H_2S$  collected is limited by the capacity of the  $AgNO_3$ -impregnated filter.
- 1.3.7. Another disadvantage is the necessity of protecting the impregnated filters from light at all times. The  $AgNO_3$ -impregnated filters darken upon exposure to light.
- 1.3.8. Alkaline cyanide solutions are used in the sample preparation and analytical procedure. Safety precautions must be followed during their use and disposal.
- 1.3.9. The alkaline cyanide solutions may contain some background sulfide, which must be subtracted from the amount of sulfide found in each sample. The precision of the method is affected by the reproducibility of the background sulfide.

Withdrawn  
Provided for Historical Reference only

1.4. Physical and chemical properties (8.4., 8.5.)

Hydrogen sulfide (CAS No. 7783-06-4) is a colorless, poisonous, flammable gas with a characteristic odor of rotten eggs at low concentrations. Collapse, coma, and death from respiratory failure may result from brief exposure at high concentrations.

Exposure at low concentrations produces irritation of conjunctiva and mucous membranes. It is soluble in water, alcohol, ether, petroleum solvents, and crude petroleum. Some physical properties are listed:

Molecular formula	-	H <sub>2</sub> S
Molecular weight	-	34.08
Specific gravity	-	1.1.92 (air = 1)
Melting point	-	-82.9 °C
Boiling point	-	61.8 °C
Auto-ignition temp.	-	250 °C
Explosive range in air	-	4.5 - 45.5%
Odor threshold	-	0.02 ppm
Olfactory fatigue level	-	100 ppm

Withdrawn  
Provided for Historical Reference only

1.5. Occupations with Potential Exposure to Hydrogen Sulfide (8.4.)

Animal fat and oil processors	Lead Removers
Animal manure removers	Lithographers
Artificial flavor makers	Lithopone
Asphalt storage workers	Livestock farmers
carbonate makers	Barium Manhole and trench workers
Blast furnace workers	Metallurgists
Brewery workers	Miners
Bromide-brine workers	Natural gas
production and Cable splicers	processing workers
Caisson workers	Painters using polysulfide
Carbon disulfide makers	caulking compounds
Cellophane makers	Papermakers
Chemical laboratory workers,	Petroleum production and
Teachers, students	refinery workers
Cistern cleaners	Phosphate purifiers
Citrus root fumigators	Photoengravers
Coal gasification workers	Pipeline maintenance
workers Coke oven workers	Pyrite burners
Copper-ore sulfidizers	Rayon makers
Depilatory makers	Refrigerant makers
Dyemakers	Rubber/plastics processors
Excavators	Septic tank cleaners
Felt makers	Sewage treatment plant
workers Fermentation process workers	Sewer workers
Fertilizer makers	Sheepdippers
Fishing and fish-processing workers	Silk makers
Fur dressers	Slaughterhouse
workers Geothermal-power drilling	Smelting workers
and production workers	Soapmakers
Gluemakers	Sugar beet and cane
processors Gold-ore workers	Sulfur spa workers
Heavy-metal precipitators	Sulfur products
processors Heavy-water manufacturers	Synthetic-fiber
workers Hydrochloric acid purifiers	Tank gagers
Hydrogen sulfide production	Tannery
workers and sales workers	Textiles printers
Landfill workers	Thiophene workers
Lead ore sulfidizers	Well diggers and
cleaners Wool pullers	

**2. Range and Detection Limit (8.6.)**

- 2.1. The analytical working range is from 0.05 to 4 µg/mL as sulfide.
- 2.2. The qualitative and quantitative detection limits of the analytical method for a 50 mL sample extraction volume are 1.0 and 2.5 µg H<sub>2</sub>S, respectively. These values correspond to 0.4 ppm and 0.90 ppm H<sub>2</sub>S when using a 2 L air volume.

### 3. Method Performance (8.6.)

- 3.1. This method was evaluated at high (85 to 88%) RH over the range of 10.1 to 39.2 ppm at an approximate temperature and atmospheric pressure of 25.5 °C and 640 mmHg, respectively. Samples were taken for 10 min using a 0.2 L/min sampling rate.
- 3.2. The pooled coefficient of variation ( $CV_T$ ) for the overall sampling and analytical method in the range tested was 0.038 and bias was -0.031.
- 3.3. In validation experiments, this method was capable of measurements within  $\pm 25\%$  of the true value at least 95% of the time. The overall error of the method was +10.7%.
- 3.4. The collection efficiency at high (85-88%) RH was 100%. A concentration of 40 ppm was used. Breakthrough tests were conducted at both high (86%) and low (18%) RH at the same concentration of 40 ppm. At either RH, 6% breakthrough of  $H_2S$  occurred after 20 min at a sampling rate of 0.2 L/min.
- 3.5. In storage stability studies, the mean recovery of samples analyzed after 30 days was within 3% of the mean recovery of similar samples analyzed immediately after collection. Sample concentrations were about 20 ppm. The samples were stored at normal lab temperatures in a dark environment.

### 4. Interferences

- 4.1. The presence of other inorganic sulfide compounds in the atmosphere will interfere with the analysis of  $H_2S$ .
- 4.2. Any substance in the atmosphere which oxidizes the sulfide formed on the impregnated filter can be a negative interference. Literature studies suggest the presence of gases such as  $SO_2$ ,  $CO_2$ ,  $O_2$ ,  $NO_2$ , and  $NH_3$  in the atmosphere will not affect the recovery of  $H_2S$  (8.2., 8.3.).
- 4.3. Any compound with the same peak potential as sulfide when using the analytical conditions described in this method is an interference. Studies suggest that organic and other inorganic sulfur compounds (sulfites, thiosulfates, mercaptans, etc.) will not interfere in the analysis (8.7., 8.8.).
- 4.4. If an interference exists, changing the operating conditions of the polarograph or the electrolyte may circumvent the problem.
- 4.5. When other substances are known or suspected to be present in the air sampled, the identities of the substances should be transmitted with the sample.

### 5. Sampling

- 5.1. Equipment
  - 5.1.1. Sample assembly:
    - Filter holder consisting of a two-piece cassette, 37-mm diameter. Backup pad, 37-mm, cellulose.
    - Silver nitrate-impregnated cellulose filters, 37-mm (see the Appendix if preparation is necessary)

Withdrawn  
Provided for Historical Reference only

- 5.1.2. Gel bands (Omega Specialty Instrument Co., Chelmsford, MA) for sealing cassettes.
- 5.1.3. Sampling pumps capable of sampling at 0.1 to 0.2 liters per minute (L/min).
- 5.1.4. Assorted flexible tubing.
- 5.1.5. Stopwatch and bubble tube or meter for pump calibration.

5.2. Sampling Procedure

**(Note: The impregnated filters should be handled with nonmetallic forceps. Also, the AgNO<sub>3</sub>-impregnated filters are light-sensitive and can turn black upon prolonged exposure to light.)**

- 5.2.1. Place an impregnated filter and a cellulose backup pad into a 37 mm polystyrene two-piece cassette filter holder. Seal the cassette with a shrinkable gel band. Place plastic end caps on the cassette. Shield the cassette from light by completely wrapping it with black tape or aluminum foil. Store the cassette in a dark environment until use.
- 5.2.2. Calibrate each personal sampling pump with a prepared cassette in-line at the flow rate listed below:

<u>Type of Sample</u>	<u>Flow Rate</u>	<u>Sampling Time</u>
Ceiling, STEL or Peak	0.2 L/min	15 min
TWA	0.1 L/min	60 min

- 5.2.3. Attach prepared cassettes to calibrated sampling pumps (the backup pad should face the pump) and place in appropriate positions on the employee or workplace area.

Begin sampling; due to potential breakthrough, **do not exceed** the recommended sampling times listed above.

- 5.2.4. Place plastic end caps on each cassette after sampling.
- 5.2.5. If possible, refrigerate the collected samples during storage periods prior to analysis. Samples do not need to be refrigerated during handling or shipping.

5.3. Shipment

Submit at least one blank sample with each set of air samples. Blank filter samples should be handled in the same manner as other samples, except no air is drawn through the blank. Attach an OSHA-21 seal around each cassette in such a way as to secure the end caps. Send the samples to the laboratory with the OSHA 91A paperwork requesting hydrogen sulfide analysis.

**6. Analysis**

6.1. Safety Precautions

- 6.1.1. Safety glasses, lab coat, and gloves must be worn at all times.

Withdrawn  
Provided for Historical Reference only

- 6.1.2. Cyanide compounds and solutions are powerful poisons which prevent the utilization of oxygen by the body tissues. Care must be exercised when using these compounds. They can be absorbed through the skin. Cyanides which contact the skin should be washed off immediately.
- 6.1.3. Mercury vapor and liquid are very toxic substances which can be readily absorbed through the respiratory tract and the skin. Avoid skin contact with mercury and clean up mercury spills immediately.
- 6.1.4. The analyst must properly dispose of all reagents, standards, prepared samples, and waste solutions after their use. Disposal of solutions containing no mercury is accomplished by pouring the solution down a sink with copious amounts of tap water and rinsing the containing vessel several times with DI H<sub>2</sub>O. Waste solutions of mercury and cyanide from the polarograph cell cups should be flushed in a large glass beaker with copious amounts of DI H<sub>2</sub>O in a sink before placing the waste mercury in a metal container specifically designed for recycling. Exercise care to prevent any mercury from entering the sink. Rinse the empty used polarograph cell cups several times with DI H<sub>2</sub>O before cleaning. Never acidify the solutions or dispose in drains which are used for acid disposal. Any acidic environment has the potential to convert the cyanide or sulfide salts to deadly hydrogen cyanide and/or H<sub>2</sub>S gas.
- 6.1.5. Waste solutions containing CdS precipitate should be poured into a labeled plastic waste bottle containing 0.1 M NaOH. The waste bottle should be properly disposed of. Check state and federal regulations for proper disposal. Never pour this waste solution down the sink.
- 6.1.6. Refer to the Standard Operating Procedure (SOP) (8.9.) and instrument manufacturer manuals for proper operation of the polarographic instrument and safety precautions.
- 6.2. Equipment
  - 6.2.1. Polarographic Analyzer or Controller, Model 384 or 384B, (Princeton Applied Research (PAR), Princeton, NJ) with a Model 303 or 303A dropping mercury electrode.
  - 6.2.2. Reference salt bridge tube (PAR Model K0154) filled with a saturated potassium chloride solution.
  - 6.2.3. Saturated calomel reference electrode (PAR Model K0077) for use with the reference electrode salt bridge tube.
  - 6.2.4. Glass or polypropylene polarographic cells, 15-mL.
  - 6.2.5. Nitrogen purification system: Gas purifier for deoxygenating nitrogen, (Oxiclear, part no. DGP-250, Labclear, Oakland, CA. As an alternative, an oxygen scrubber can be constructed using a vanadous chloride solution as described in reference 8.10.).
  - 6.2.6. Sulfide specific ion electrode (Model 94-16A, Orion Research, Inc. Cambridge, MA).
  - 6.2.7. Double junction reference electrode (Model 90-02, Orion Research) with inner and



Withdrawn  
Provided for Historical Reference only

outer filling solutions (Cat. Nos. 90-00-02 and 90-00-03, respectively).

- 6.2.8. Millivolt meter.
  - 6.2.9. Magnetic stirrer and Teflon stirring bars.
  - 6.2.10. Burette, glass, 10-mL.
  - 6.2.11. Beakers, polypropylene, 100-mL, with tight fitting covers.
  - 6.2.12. Bottle, narrow mouth, amber, linear polyethylene, opaque, 125- or 250-mL, with screw cap - used for storage of light-sensitive material.
  - 6.2.13. Bottles, wide mouth, polypropylene, with screwcaps.
  - 6.2.14. Nonmetallic forceps.
  - 6.2.15. Glass rod.
  - 6.2.16. Polystyrene disposable beakers.
  - 6.2.17. Automatic pipette, adjustable, 20- to 5000- $\mu$ L range (Gilson Pipetman, Rainin Instruments Inc.).
  - 6.2.18. Volumetric pipettes, volumetric flasks, beakers and other laboratory glassware.
  - 6.2.19. Tissue paper.
  - 6.2.20. Analytical balance (0.01 mg).
- 6.3. Reagents: All chemicals should be reagent grade or better.
- 6.3.1. Deionized water (DI H<sub>2</sub>O): Deionized, filtered, deoxygenated water is needed for preparation of all solutions, which will be used in the analysis. It is recommended to purge the DI H<sub>2</sub>O with purified nitrogen before use.
  - 6.3.2. Sodium hydroxide (NaOH).
  - 6.3.3. Sodium hydroxide, 0.1 M solution. Dissolve 4.0 g of N in approximately 400 mL of DI H<sub>2</sub>O and then dilute the solution to 1 L. Store the solution in a wide-mouth polypropylene bottle.
  - 6.3.4. Sodium cyanide (NaCN), containing s 0.001% sulfide as an impurity.
  - 6.3.5. Alkaline 0.5 M NaCN solution. Dissolve 2.0 g of NaOH and then 12.25 g of NaCN in approximately 300 mL of DI H<sub>2</sub>O and then dilute the solution to 500 mL. Immediately and carefully transfer the solution to a wide-mouth, polypropylene bottle.
  - 6.3.6. Alkaline 0.25 M NaCN solution. Dilute 250 mL of the 0.5 M NaCN solution to 500 mL with 0.1 M NaOH solution. Immediately and carefully transfer the solution to a wide-mouth polypropylene bottle.
  - 6.3.7. Nitric acid (HNO<sub>3</sub>), concentrated.
  - 6.3.8. Nitric acid, 6 M: Carefully dilute 384 mL of concentrated HNO<sub>3</sub> to 1 L using DI H<sub>2</sub>O.
  - 6.3.9. Nitric acid, 1 M: Dilute 6.4 mL of concentrated HNO<sub>3</sub> to 100 mL with DI H<sub>2</sub>O.
  - 6.3.10. Sodium sulfide nonahydrate (Na<sub>2</sub>S·9H<sub>2</sub>O), crystalline.

Withdrawn  
Provided for Historical Reference only

6.3.11. Sodium sulfide ( $\text{Na}_2\text{S}$ ) stock solution (equivalent to 1,000  $\mu\text{g/mL}$   $\text{H}_2\text{S}$ ):

1. Weigh 0.7050 g of  $\text{Na}_2\text{S}\cdot 9\text{H}_2\text{O}$  crystals in a 50 mL disposable polystyrene beaker (Note: Before weighing, remove any surface impurities by rinsing the crystals with DI  $\text{H}_2\text{O}$  in an exhaust hood. Quickly blot the crystals dry with tissue paper. Handle the crystals with nonmetallic forceps.).
2. Dissolve the weighed  $\text{Na}_2\text{S}\cdot 9\text{H}_2\text{O}$  crystals in a disposable beaker containing about 30 mL of 0.1 M NaOH.
3. Quantitatively transfer the solution to a 100 mL volumetric flask. Rinse the beaker several times with 0.1 M NaOH and transfer to the volumetric flask. Dilute the flask to the mark with 0.1 M NaOH.
4. Wrap the flask with aluminum foil to protect it from light. Refrigerate and keep the flask in a dark environment prior to and immediately after each use.
5. Prepare this solution weekly.

6.3.12. Cadmium sulfate ( $\text{CdSO}_4$ ) anhydrous.

6.3.13. Cadmium sulfate ( $\text{CdSO}_4$ ), 0.1500 M standard solution (used in the standardization of the  $\text{Na}_2\text{S}$  stock solution).

1. Heat approximately 3.5 g of anhydrous  $\text{CdSO}_4$  at 110 °C for 24 h in a drying oven. Cool the dried  $\text{CdSO}_4$  in a desiccator for approximately 30 min.
2. Quickly weigh 3.127 g of the dried  $\text{CdSO}_4$  in a 50 mL disposable polystyrene beaker and then dissolve in approximately 25 mL of DI  $\text{H}_2\text{O}$ .
3. Quantitatively transfer this solution to a 100 mL volumetric flask. Rinse the beaker several times with DI  $\text{H}_2\text{O}$  and transfer the rinses to the volumetric flask. Dilute to volume with DI  $\text{H}_2\text{O}$ .

6.3.14. Saturated potassium chloride solution for filling the salt bridge.

6.4. Standard Preparation

6.4.1. To standardize the  $\text{Na}_2\text{S}$  stock solution a titration is performed with a standard 0.1500 M  $\text{CdSO}_4$  solution. A sulfide ISE is used as the indicator. An orange-colored precipitate ( $\text{CdS}$ ) is formed from the reaction of  $\text{Cd}^{2+}$  with  $\text{S}^{2-}$  (8.11.).

1. Connect the sulfide ion specific and double junction reference electrodes to the millivolt meter. Place the electrodes in a 0.1 M NaOH solution to reduce stabilization time.
2. Allow the  $\text{Na}_2\text{S}$  stock solution to come to ambient temperature and then pipette 25 mL of this solution into a 50 mL disposable beaker containing a Teflon stirring bar. Use magnetic stirring throughout the titration.
3. Place the electrodes in the beaker and record the initial mV reading. Place a piece of insulating material such as cardboard or Styrofoam between the magnetic stirrer and the beaker to minimize the heating of the standard solution by the stirrer.
4. Record the electrode potential as a function of titrant volume added. Initially add the 0.1500 m  $\text{CdSO}_4$  standard titration solution from a 10 mL burette in 0.5 to 1.0 mL increments to the beaker. When the potential change per increment begins to

Withdrawn  
Provided for Historical Reference only

increase, add 0.05 mL to 0.2 mL of the titrant to the beaker. Add about 1 mL of the titrant beyond the endpoint.

5. The end point is the point of greatest slope on the titration curve. This can be determined by plotting the electrode potential versus titrant volume or tabulating the change in potential per mL increment of titrant added.
6. Calculate the  $\mu\text{g/mL}$  equivalent of  $\text{H}_2\text{S}$  in the  $\text{Na}_2\text{S}$  stock solution:

$$\text{H}_2\text{S } (\mu\text{g/mL}) = \frac{A \times B \times 32.064 \times 1,000 \times 1.0629}{C}$$

Where:

- A = molarity of the standard  $\text{CdSO}_4$  solution.
- B = mL of standard  $\text{CdSO}_4$  solution required to titrate the  $\text{Na}_2\text{S}$  stock solution.
- 32.064 = milligrams per milliequivalent sulfide.
- 1,000 = micrograms per milligram
- 1.0629 = gravimetric factor ( $\text{H}_2\text{S}/\text{S}^{2-}$ )
- C = 25.0 mL = mL of the  $\text{Na}_2\text{S}$  solution used.

- 6.4.2. Prepare standards equivalent to 100 and 10  $\mu\text{g/mL}$   $\text{H}_2\text{S}$  by making appropriate serial dilutions of the  $\text{Na}_2\text{S}$  stock solution with the 0.1 M NaOH solution. Protect these standards from light and prepare daily.

#### 6.5. Sample Preparation

Wash all glassware and plastic ware, excluding the polarograph cell cups, in detergent, rinse in dilute (1 M) nitric acid, rinse thoroughly with DI  $\text{H}_2\text{O}$ , and then air dry prior to use. Wash the cell cups in 6 M  $\text{HNO}_3$  and rinse with DI  $\text{H}_2\text{O}$  before use.

- 6.5.1. Open the filter cassette. Carefully remove the sample filter with nonmetallic forceps and place it in a 100 mL polypropylene beaker.
- 6.5.2. Immediately pipette 25 mL of the alkaline 0.5 M NaCN solution into the beaker and then tightly cover it with a polypropylene lid.
- 6.5.3. Place the covered beaker in a dark environment at ambient room temperature for approximately 6 h. Occasionally swirl the sample solution with a gentle motion. The silver sulfide on the filter will dissolve in the cyanide solution during this time.
- 6.5.4. Store the covered sample beaker overnight in a refrigerator.
- 6.5.5. Remove the beaker from the refrigerator and place in a dark environment for approximately 45 min to allow the sample solution to equilibrate to room temperature.
- 6.5.6. Pipette 25 mL of the 0.1 M NaOH solution into the sample beaker and then stir the sample solution with a clean glass rod. Cover each sample solution with a plastic lid until analysis. The solutions should be analyzed sometime during the same work day.

Withdrawn  
Provided for Historical Reference only

6.6. Analytical Procedure

- 6.6.1. Protect and electrically isolate the silver reference electrode wire on the DME by installing a clean, empty reference electrode glass sleeve with a Vycor tip.
- 6.6.2. Install the filled salt bridge tube on the electrode support block. Insert the saturated calomel reference electrode (SCE) into the salt bridge tube and then connect the SCE cable to the reference electrode input.
- 6.6.3. Refer to reference 8.9. for operating the polarograph.
- 6.6.4. Turn on the nitrogen to a pressure of about 2 psi.
- 6.6.5. Initially set the following values (Note: Some of the settings mentioned below are instrument specific. Refer to specific operating and service manuals for other types of polarograph):

Analytical Technique	-	DPP
Initial Potential	-	-0.610 v
Final Potential	-	-0.770 v
Peak Location	-	yes
Peak Potential	-	-0.690 V
Purge Time	-	60 s
Drop Time	-	1 s
Scan Increment	-	2 mV Scan
Rate	-	2mV/s
Pulse Height	-	0.050 V
Tangent Fit	-	No Derivative
No Blank Subtraction	-	Yes
Override Position	-	10 11 12 13 14 15 16 17
Status	-	Y N Y Y N Y Y N
Drop Size	-	Large

- 6.6.6. The differential pulse polarogram of sulfide using a DME yields a peak at approximately -0.690 V. The peak potential depends on reference electrode response and the DME operating parameters.
- 6.6.7. Prepare standards in the range of 0.05 to 4.0  $\mu\text{g/mL}$  immediately before analysis by diluting  $\mu\text{L}$  aliquots of the  $\text{Na}_2\text{S}$  working standards to 10 ml in the polarograph cell cups with the alkaline 0.25 M NaCN solution. A suggested method of preparation of the standard solutions in the cell cups is given:

Withdrawn  
 Provided for Historical Reference only

SS H <sub>2</sub> S Concn* (µg/mL)	Aliquot of SS (µL)	0.25 M NaCN (ML added)	Final H <sub>2</sub> S Concn (µg/mL)
0**	---	9.80	Reagent Blank
10	50	9.95	0.050
10	100	9.90	0.100
10	200	9.80	0.200
100	50	9.95	0.500
100	100	9.90	1.00
100	200	9.80	2.00
100	300	9.70	3.00
100	400	9.60	4.00

\* All working standards are prepared in 0.25 M NaCN/0.1 N NaOH solution.

\*\* Reagent Blank contains 9.8 mL NaCN/0.1 N NaOH and 0.2 mL of 0.1 N NaOH. SS = Na<sub>2</sub>S stock solution. All H<sub>2</sub>S concns listed are as equivalent concns.

- 6.6.8. Scan a 1 µg/mL standard to determine the actual sulfide peak potential. If the peak is significantly different from -0.69 V, reset the peak, initial, and final potentials to window the peak (±0.080 V from peak).
- 6.6.9. Observe in this polarogram the rapidly rising baseline when at a more positive potential than the sulfide peak. This rising baseline is due to cyanide in the supporting electrolyte, which has a peak potential of approximately -0.36 V. This baseline is minimized by a reagent blank subtraction routine found in the polarograph software.
- 6.6.10. Pipette 10 mL aliquots of the 50 mL sample solutions from the polypropylene beakers into the polarograph cell cups immediately prior to analysis.
- 6.6.11. Record the peak potential (V) and the peak current (µA) for each standard and sample.
- 6.7. Analytical Recommendations
  - 6.7.1. Verify that mercury is flowing properly through the glass capillary of the DME before the analysis begins and during the analysis by observing mercury drop formation.
  - 6.7.2. Because the sensitivity and linearity range increase with increasing cyanide concentration, it is important the cyanide concentration in the standard and analytical solutions be identical.
  - 6.7.3. Analyze the reagent blank, a series of standards in the range of interest, and the samples. A freshly prepared standard should be analyzed after every five or six samples to monitor polarograph performance. Always rinse the electrodes with DI H<sub>2</sub>O and blot dry with a clean dry tissue before analyzing the next solution.

Withdrawn  
Provided for Historical Reference only

## 7. Calculations

7.1. Determine the  $\mu\text{g/mL}$   $\text{H}_2\text{S}$  content of each sample and blank using a concentration-response linear regression curve ( $\mu\text{A}$  vs.  $\mu\text{g/mL}$ ).

7.2. Calculate the total micrograms contained in each sample:

$$\mu\text{g H}_2\text{S} = [(\text{Sample Vol, mL}) \times \mu\text{g/mL H}_2\text{S}]$$

7.3. Each air sample is blank corrected ( $\mu\text{g}$  sample -  $\mu\text{g}$  blank) and the concentration is then calculated to determine  $\text{H}_2\text{S}$  exposure using the following equation:

$$\text{ppm H}_2\text{S} = \frac{\text{MV} \times \mu\text{g H}_2\text{S}_b}{\text{Molecular Weight} \times (\text{Air Volume, L})}$$

Where:

$$\text{MV (Molar Volume)} = 24.45 \text{ (25 }^\circ\text{C and 760 mmHg)}$$

$$\mu\text{g H}_2\text{S}_b = \text{Blank corrected sample result}$$

$$\text{Molecular Weight (H}_2\text{S)} = 34.080$$

This equation reduces to:

$$\text{ppm H}_2\text{S} = \frac{0.717 \times \mu\text{g Sample}}{(\text{Air Volume, L})}$$

7.4. Reporting Results

Results are reported to the industrial hygienist as ppm  $\text{H}_2\text{S}$ .

## 8. References

8.1. **National Institute for Occupational Safety and Health:** NIOSH Manual of Analytical Methods, 2nd ed., Vol. 2 (DHEW/NIOSH Pub. No. 77-157-B). Cincinnati, OH: National Institute for Occupational Safety and Health, 1977. pp. S4- 1-S4-10.

8.2. **Natusch, D.F.S., H.B. Klonis, H.D. Axelrod, R.D. Teck and J.P. Lodge, Jr.:** Sensitive Method for Measurement of Atmospheric Hydrogen Sulfide. Anal. Chem. **44:** 2067-2070 (1972).

8.3. **Natusch, D.F.S., J.R. Sewell and R.L. Tanner:** Determination of Hydrogen Sulfide in Air -- An Assessment of Impregnated Paper Tape Methods. Anal. Chem. **46:** 410- 415 (1974).

8.4. **National Institute for Occupational Safety and Health:** Criteria for a Recommended Standard - Occupational Exposure to Hydrogen Sulfide (DHEW/NIOSH Pub. No. 77-158). Cincinnati, OH: National Institute for Occupational Safety and Health, 1977.

8.5. **Fassett, D.W. and D.D. Irish, ed.:** Patty's Industrial Hygiene and Toxicology, 2nd rev. ed., Vol. 2. New York: John Wiley and Sons, 1963.

8.6. **Occupational Safety and Health Administration Technical Center:** Hydrogen Sulfide Backup Data Report (ID-141). Salt Lake City, UT. Revised 1989.

8.7. **Renard, J.J., G. Kubes and H.I. Bolker:** Polarographic Determination of Sulfur Compounds in Pulping Liquors. Anal. Chem. **47:** 1347-1352 (1975).

Withdrawn  
Provided for Historical Reference only

- 8.8. **Noel, D.L.:** Sulfur Compounds in Kraft Pulping Liquor. Tappi 61: 73-76 (1978).
- 8.9. **Occupational Safety and Health Administration Technical Center:** Standard Operating Procedure for Polarography. Salt Lake City, UT. In progress (unpublished).
- 8.10. **Princeton Applied Research:** Application note 108, Why Dearthation... and How. Princeton, NJ: Princeton Applied Research, 1974.
- 8.11. **Ehman, D.L.:** Determination of Parts-per-Billion Levels of Hydrogen Sulfide in Air by Potentiometric Titration with a Sulfide Ion-Selective Electrode as an Indicator. Anal. Chem. 48: 918-920 (1976).

**Appendix**

Preparation of H<sub>2</sub>S Sampling Devices - Impregnated Filters

---

Preparation of filter impregnation solution, 2% AgNO<sub>3</sub>

---

Reagents (All reagents should be reagent grade or better):

Glycerin.

Ethanol, absolute.

Silver nitrate (AgNO<sub>3</sub>).

---

1. Dissolve 2 g AgNO<sub>3</sub> in approximately 60 mL DI H<sub>2</sub>O contained in a glass beaker. Immediately transfer this solution to a 100-mL volumetric flask.
  2. Pipette 1 mL of 1 N HNO<sub>3</sub>, 5 mL of glycerin, and 20 mL of ethanol into the 100-mL volumetric flask. Dilute the solution to volume with DI H<sub>2</sub>O.
  3. Immediately transfer the AgNO<sub>3</sub> solution to an amber polyethylene bottle. Store the bottle in a dark environment. The AgNO<sub>3</sub> solution has at least a 6 month shelf-life.
- 

Preparation of Impregnated Filters

---

Whatman No. 4 (Special Order, Whatman Labsales, Hillsboro, OR) cellulose filters (37-mm diameter) are impregnated with the AgNO<sub>3</sub> solution described above. This filter is prepared as follows:

Withdrawn  
Provided for Historical Reference only

---

---

1. Place the cellulose filter over one of the holes of a polypropylene test tube rack or scintillation vial support rack. The rack should have hole diameters of 30-mm.
  2. Using subdued light, pipette 250  $\mu\text{L}$  of the  $\text{AgNO}_3$  impregnation solution on the filter.
  3. Carefully transfer the filters and support to a dark environment and allow to dry at ambient temperature for approximately 120 min.
  4. Handle the filters in subdued light.
  5. Store the filters together in a dark environment at normal lab temperatures in tightly covered polystyrene petri dishes. Filters can be stored for at least 3 months.
- 
-



## Hydrogen Sulfide

OSHA Method ID-141 Backup

Revised December 1989

### Introduction

The procedure for collection and analysis of air samples for hydrogen sulfide (H<sub>2</sub>S) is described in OSHA Method No. ID-141 (12.1). Briefly, H<sub>2</sub>S is collected on Whatman No. 4 filter paper (Special Order, Whatman Labsales, Hillsboro, OR) which has been impregnated with silver nitrate (AgNO<sub>3</sub>). The H<sub>2</sub>S reacts with AgNO<sub>3</sub> to form silver sulfide, a greyish-black precipitate (12.2, 12.3). The silver sulfide is dissolved in an alkaline cyanide solution and then analyzed by differential pulse polarography (DPP) for sulfide.

### 1. Experimental Procedure

This method has been evaluated using a sampling rate of 0.2 L/min and a total air volume of 2 L. The OSHA Ceiling Permissible Exposure Limit (PEL) was 20 ppm when this evaluation was performed; therefore, concentrations of approximately 10, 20, and 40 ppm H<sub>2</sub>S were used for the evaluation. The evaluation followed the NIOSH protocol (12.4) for method evaluation with some exceptions and consisted of the following experiments:

1. An analysis of a total of 18 spiked samples (6 samples at each of three test levels) prepared by spiking known amounts of sodium sulfide.
2. An analysis of a total of 18 samples (6 samples at each of the three test levels) collected from dynamically generated test atmospheres using a 2-L total air volume.
3. A determination of the collection efficiency and breakthrough of the AgNO<sub>3</sub>-impregnated filters at about 40 ppm.
4. A test of the storage stability of 18 collected samples.
5. A determination of the detection limit of the method.
6. A comparison of sampling and analytical methods.
7. An assessment of the precision and accuracy of the method.
8. An assessment of the method using an independent laboratory.
9. Conclusions - including a discussion of recent PEL changes for H<sub>2</sub>S.

---

Note: During this evaluation, all AgNO<sub>3</sub>-impregnated filter samples were taken using a sample collection flow rate of approximately 0.2 L/min. All generated samples were taken at an atmospheric pressure of approximately 640 mmHg and a temperature of 25.5 °C with the exception of the independent laboratory test (Section 10). These samples were taken at atmospheric conditions of approximately 761 mmHg and 21 °C.

---

### 2. Analysis

Samples were spiked and analyzed to determine analytical performance of the method and to assess desorption efficiency.

# Withdrawn

## Provided for Historical Reference only

**Procedure:** Samples were prepared by adding known amounts of a standardized sodium sulfide ( $\text{Na}_2\text{S}$ ) stock solution to  $\text{AgNO}_3$ -impregnated Whatman 4 filters. Approximately 25, 50, and 100  $\mu\text{g}$  (as  $\text{H}_2\text{S}$ ) were spiked. Details are discussed below.

- 2.1 A  $\text{Na}_2\text{S}$  stock solution was prepared by dissolving and diluting 35.855 g of  $\text{Na}_2\text{S} \cdot 9\text{H}_2\text{O}$  to 500 mL with deoxygenated 0.1 M  $\text{NaOH}$ . The stock was standardized using the procedure listed in the method (12.1). This stock solution was equivalent to 10,175  $\mu\text{g}/\text{mL}$   $\text{H}_2\text{S}$ . The solution was protected from light and stored in a refrigerator.
- 2.2 Silver nitrate-impregnated Whatman 4 filters were prepared as mentioned in reference 12.1.
- 2.3 The three sets of spiked samples were prepared by injecting known volumes of the standardized  $\text{Na}_2\text{S}$  stock solution on the  $\text{AgNO}_3$ -impregnated filters. The spiked filters were air dried in a dark environment for approximately 24 h prior to resorption. Each set contained six samples and a sample blank.
- 2.4 Sample solutions were prepared and analyzed for sulfide by DPP as described in reference 12.1 with one slight modification: Analytical standards were prepared by serial dilution of the standardized  $\text{Na}_2\text{S}$  stock solution (10,175  $\mu\text{g}/\text{mL}$   $\text{H}_2\text{S}$ ) prepared in Section 2.1. Each sample was analyzed twice.

**Results:** The results of the spiked filter samples are given in Table 1. The desorption efficiency (DE) for each sample set analyzed was 0.923, 0.892 and 0.882, respectively. The DE at each level is significantly different from 1.0. The slightly low DE appeared to be due to difficulties in spiking the filters with the  $\text{Na}_2\text{S}$  stock solution and not to problems in the desorption of analyte. Potential problems when spiking these solutions were noted:

1. Sodium sulfide is an unstable compound (12.5). Exposure of sodium sulfide crystals to air produces  $\text{H}_2\text{S}$ .
2. Aqueous solutions of  $\text{Na}_2\text{S}$  are slowly converted to sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3$ ) from standing in contact with air.

### 3. Sampling and Analysis

**Procedure:** To determine the ability of the method to sample and analyze for  $\text{H}_2\text{S}$ , known concentrations of generated samples were prepared. Test atmospheres of  $\text{H}_2\text{S}$  gas in this and all following sections were dynamically generated by diluting  $\text{H}_2\text{S}$  gas from a cylinder with purified humid air in a Teflon mixing tee. The cylinder concentration had been certified as  $220 \pm 5$  ppm  $\text{H}_2\text{S}$  in nitrogen (Airco, Murray Hill, NJ). To verify the concentration of the cylinder, samples were taken and analyzed using the NIOSH sampling and analytical method no. S4 (12.6) with some modifications. The procedure followed is listed below.

- 3.1 Three samples of  $\text{H}_2\text{S}$  gas from the cylinder were collected in midjet impingers containing 10 mL of a cadmium hydroxide/STRactan absorbing solution at approximately 0.1 L/min for 5 min. A backup impinger containing this solution was connected in series to each sample to assure no loss of analyte from breakthrough.
- 3.2 The samples were analyzed for sulfide content by the methylene blue calorimetric procedure (12.6). Standards were prepared from the  $\text{Na}_2\text{S}$  stock solution prepared in Section 2.
- 3.3 No  $\text{H}_2\text{S}$  was detected in the backup impinger of each sample (detection limit = 1.6 ppm  $\text{H}_2\text{S}$ ).
- 3.4 The average  $\text{H}_2\text{S}$  concentration of the gas cylinder was  $208 \pm 14$  ppm (95% confidence level). This experimental concentration did not significantly differ from the manufacturers' stated

# Withdrawn

## Provided for Historical Reference only

concentration of  $220 \pm 5$  ppm. The stated value was used when calculating all theoretical (taken) results.

Test atmospheres were then generated using this gas and the generation system described in the Appendix. Six samples were simultaneously collected from the manifold for approximately 10 min. This was performed for each of the three test levels. The samples were collected on the  $\text{AgNO}_3$ -impregnated filters with model C-210 portable sampling pumps (Mine Safety Appliances, Pittsburgh, PA) at high (85 to 88%) RH and  $25.5^\circ\text{C}$ .

**Results:** The three sets of samples were analyzed twice by DPP using procedures described in the method (12.1). The results of the sampling and analysis experiment are shown in Table 2. A correction for DE was not performed on the results. The theoretical (taken) concentration of the generated gas was determined from the measured flow rates of the diluting air and the certified  $\text{H}_2\text{S}$  gas from the cylinder at each test level. Two results were deleted as outliers since they failed to pass the outlier test at the 99% confidence level.

#### 4. Collection Efficiency (CE)

**Procedure:** Six samples were simultaneously collected on  $\text{AgNO}_3$ -impregnated filters which were connected in series to backup midget impingers containing 10 ml of a cadmium hydroxide/STRactan absorbing solution. Samples were collected at high (85-88%) RH and  $25.4^\circ\text{C}$ . The samples were collected at a concentration of approximately 40 ppm for 10 min [Note: These were the same samples which were generated and analyzed in the sampling and analysis experiment (Section 3) at about 40 ppm]. The amounts of  $\text{H}_2\text{S}$  collected in the filters and backup impingers were then measured.

**Results:** The CE of each filter sample was calculated by dividing the amount of  $\text{H}_2\text{S}$  collected on the filter by the total amount of  $\text{H}_2\text{S}$  collected in the filter and backup impinger. The results are given in Table 3. The CE was 100%.

#### 5. Breakthrough

**Procedure:** Two samples were simultaneously collected at approximately 40 ppm  $\text{H}_2\text{S}$  for 10, 15, and 20 minutes. The generation system test atmospheres were produced at a low (18%) and then a high (86%) RH ( $25^\circ\text{C}$ ) to determine any humidity effect on breakthrough. Each sample was collected on a  $\text{AgNO}_3$ -impregnated filter connected in series to a backup midget impinger containing 10 ml of cadmium hydroxide/STRactan absorbing solution. The midget impinger samples were analyzed for sulfide content using NIOSH method no. S4 (12.6) for  $\text{H}_2\text{S}$  with some modifications. Breakthrough of  $\text{H}_2\text{S}$  into the impinger solutions at a concentration greater than 5% of the amount generated was considered significant.

**Results:** Breakthrough was determined for each sampling period by dividing the average amount of  $\text{H}_2\text{S}$  collected in the backup impingers by the known generated  $\text{H}_2\text{S}$  concentration. The results are presented in Table 4. Breakthrough of 6% occurred at both humidities and for a sampling period of 20 min. Breakthrough was below 5% for shorter sampling periods.

#### 6. Storage Stability

A study was conducted to assess the stability of  $\text{H}_2\text{S}$  air samples collected on the  $\text{AgNO}_3$ -impregnated filters.

**Procedure:** Three sets (6 samples and a blank sample in each set) of  $\text{H}_2\text{S}$  samples were generated using  $\text{AgNO}_3$ -impregnated filters as the collection media. The generation system was set at high (85-88%) RH and  $25.5^\circ\text{C}$ . The samples were collected at approximately 20 ppm for 10 min. The samples were then stored in a dark environment at normal laboratory temperatures. A set of samples was analyzed after sample storage of 5, 14, and 30 days.

**Results:** The results of the storage stability study are shown in Table 5 [Note: The six samples listed for the 0 day storage test are the same samples taken at 20 ppm for the sampling and analysis (Section 3) experiment]. The results indicate samples are stable for a period of at least 30 days when stored under normal lab temperatures and in a dark environment.

## 7. Detection Limit

Procedure: Standard solutions in 0.1 N NaOH were prepared by serial dilutions of a stock standard solution of Na<sub>2</sub>S. The procedure used is identical to the preparation of working standards in the method (12.1), with the exception that the concentrations used for this experiment were from 10.2 to 203.5 ng/mL. Six standards at each concentration and six reagent blanks were prepared and analyzed.

### Results:

#### 7.1 Qualitative Detection Limit

The analytical procedure requires the use of a blank subtraction software routine for each sample or standard. This routine sets any constant background signal to zero. Therefore, parametric or non-parametric tests which require measurable analytical signals from reagent blank samples could not be used to determine the qualitative detection limit (12.7). The qualitative limit was estimated from results in Table 6 to be 0.020 µg/mL or 1.0 µg of H<sub>2</sub>S in a 50 ml sample extraction volume. This corresponds to 0.4 ppm H<sub>2</sub>S (2-L air volume). The coefficient of variation (CV) of replicate determinations of standards at this level was 0.128. Concentrations of H<sub>2</sub>S below 0.020 µg/mL were not consistently resolved when analyzed.

#### 7.2 Quantitative Detection Limit

The quantitative detection limit is 0.050 µg/mL or 2.5 µg of H<sub>2</sub>S in a 50 ml sample extraction volume. This corresponds to 0.9 ppm H<sub>2</sub>S (2-L air volume). As shown in Table 6, the CV of replicate determinations of standards at this concentration is less than 0.10.

## 8. Method Comparison

A side-by-side method comparison was performed as an independent measurement of the generated concentrations. The NIOSH sampling and analytical method no. S4 (12.6) was used with some modifications. The samples for the NIOSH method were taken simultaneously with impregnated filter samples during the sampling and analysis experiment (Section 3).

**Procedure:** Six samples (NIOSH method) of the test atmosphere were simultaneously collected from the sampling manifold for approximately 10 min at a flow rate of 0.2 L/min at each of the three test levels. Each set of six samples was concurrently collected with the filter samples at that test concentration. The procedures used for the NIOSH method are listed below.

8.1 Samples were collected in midget impingers containing 10 ml of a cadmium hydroxide/STRactan absorbing solution with MSA Model C-210 portable sampling pumps.

8.2 The samples were analyzed for their sulfide content by the methylene blue colorimetric procedure with standards prepared from the standardized Na<sub>2</sub>S stock solution.

**Results:** The average H<sub>2</sub>S concentration (95% confidence level) determined in each set of impinger samples is:

Withdrawn  
Provided for Historical Reference only

H <sub>2</sub> S Taken Concentration	H <sub>2</sub> S Found Concentration
10.06 ppm	10.07 ± 0.17 ppm
22.18 ppm	21.94 ± 0.29 ppm
39.22 ppm	37.19 ± 0.92 ppm

The first two experimental concentrations found did not significantly differ from the taken concentrations. The difference between the taken and found concentrations for the third test (39.22 vs. 37.19 ppm) was possibly due to analytical difficulties in recovering the cadmium sulfide precipitate from the impingers at this test level.

#### 9. Precision and Accuracy

The precision and accuracy data, based on the NIOSH statistical protocol (12.6), are presented in Tables 1 and 2. The pooled coefficients of variation for spiked (CV1) and generated (CV2) samples and the overall pooled CVT are as follows:

$$CV1 = 0.031 \quad CV2 = 0.036 \quad CVT = 0.038$$

The bias was -3.1% and overall error was ±10.7%.

#### 10. Independent Assessment of Method

**Procedure:** A series of sample filters impregnated with AgNO<sub>3</sub> were prepared at the OSHA Analytical Laboratory and submitted to an independent laboratory for sample collection at low and high humidity using their H<sub>2</sub>S generation system. Their system dynamically generated H<sub>2</sub>S gas by diluting a gas stream from a cylinder of H<sub>2</sub>S with a stream of air from a compressed air cylinder. The dilution took place in a Teflon manifold. Flow rates were controlled using rotameters. The high humidity experiment was conducted by bubbling the diluting air in water before mixing. The known (taken) concentrations were determined by measuring aliquots of the generated atmospheres with a model 5700 gas chromatography (Hewlett-Packard, Avondale, PA) equipped with a flame photometric detector. The detector response was calibrated using a H<sub>2</sub>S permeation tube.

Samples were taken at three different concentrations and two different humidities by employees of the independent laboratory. A flow rate of 0.2 L/min was used with model 222-3 sampling pumps (SKC Inc., Eighty Four, PA). The samples were then submitted to the OSHA Analytical Laboratory for analysis. The samples taken at low humidity were collected and analyzed within 30 days of preparation; high humidity samples were collected about 45 days after preparation. The samples taken at high humidity were analyzed about 70 days after the filters were impregnated.

**Results:** Results are shown in Table 7. The results of the two different humidity tests indicate good agreement with the theoretical values with the exception of samples collected at about 40 ppm. All recoveries for the 40 ppm samples were about 50% lower than expected. The independent laboratory indicated monitoring problems occurred when performing the first test (low humidity) at this concentration. The problem was attributed to a faulty H<sub>2</sub>S permeation tube when calibrating the gas chromatography. For the second high concentration test (high humidity), recoveries were also about 50% low; however, the independent laboratory indicated the impregnated filters used during the 40 ppm test were not handled according to specifications. Blank air was collected through these filters a month prior to their use.

The length of time between preparation and analysis for high humidity samples (at concentrations at or below 30 ppm H<sub>2</sub>S) indicates a storage stability of at least 45 days.

Withdrawn  
Provided for Historical Reference only

## 11. Conclusions

This sampling and analytical method has been shown to be precise and accurate for determining Ceiling exposures of 10 to 40 ppm when using 0.2 L/min flow rates for 10 to 15 min. To determine compliance with the Final Rule STEL of 15 ppm H<sub>2</sub>S, the same sampling and analytical conditions can be used.

Storage stability did not pose a significant problem under the conditions tested. Breakthrough was evident after 20 min of sampling at 40 ppm which places a limitation on sampling. Due to the potential for breakthrough, it is recommended to sample for TWA exposures at a lower flow rate of 0.1 L/min. Eight 1-h samples are recommended for TWA assessments of H<sub>2</sub>S. Samples taken near the TWA PEL of 10 ppm (0.1 L/min for 1 h) will have the same total amount collected as the samples collected at 20 ppm (0.2 L/min for 15 min). Therefore, at the lower flow rate the method should not exhibit any significant sample collection problems when determining TWA exposures.

## 12. References

- 12.1 Occupational Safety and Health Administration Technical Center: Hydrogen Sulfide in Workplace Atmospheres by T. Wilczek (USDOL/OSHA-SLTC Method No. ID-141). Salt Lake City, UT. Revised 1989.
- 12.2 Natusch, D.F.S., H.B. Klonis, H.D. Axelrod, R.D. Teck, and J.P. Lodge, Jr.: Sensitive Method for Measurement of Atmospheric Hydrogen Sulfide. *Anal. Chem.* 44: 2067-2070 (1972).
- 12.3 Natusch, D. F. S., J.R. Sewell, and R.L. Tanner: Determination of Hydrogen Sulfide in Air -- An Assessment of Impregnated Paper Tape Methods. *Anal. Chem.* 46: 410-415 (1974).
- 12.4 National Institute for Occupational Safety and Health: Documentation of the NIOSH Validation Tests by D. Taylor, R. Kupel and J. Bryant (DHEW/NIOSH Pub. No. 77-185). Cincinnati, OH: National Institute for Occupational Safety and Health, 1977.
- 12.5 Windholz, H., ed.: *The Merck Index*. 9th ed. Rahway, NJ: Merck & Co., Inc., 1976.
- 12.6 National Institute for Occupational Safety and Health: NIOSH Manual of Analytical Methods. 2nd ed., Vol. 2 (DHEW/NIOSH Pub. No. 77-157-B). Cincinnati, OH: National Institute for Occupational Safety and Health, 1977. pp. S4-1-S4-100.
- 12.7 Occupational Safety and Health Administration Analytical Laboratory: OSHA Analytical Methods Manual (USDOL/OSHA-SLTC Method No. ID-102 (Backup Report)). Cincinnati, OH: American Conference of Governmental Industrial Hygienists (Pub. No. ISBN: 0-936712-66-X), 1985.

Withdrawn  
 Provided for Historical Reference only

Table 1  
 Analysis - Hydrogen Sulfide

-----0.5 × PEL*-----			-----1 × PEL*-----			-----2 × PEL*-----		
µg Taken	µg Found	DE	µg Taken	µg Found	DE	µg Taken	µg Found	DE
25.44	22.77	0.895	50.88	43.37	0.852	101.8	85.57	0.841
25.44	23.28	0.915	50.88	44.89	0.882	101.8	88.53	0.870
25.44	22.82	0.897	50.88	46.71	0.918	101.8	89.96	0.884
25.44	22.92	0.901	50.88	46.26	0.909	101.8	91.92	0.903
25.44	24.27	0.954	50.88	45.59	0.896	101.8	92.74	0.911
25.44	24.89	0.978	50.88	45.58	0.896	101.8	90.06	0.885
N		6			6			6
Mean		0.923			0.892			0.882
Std Dev		0.035			0.023			0.030
CV1		0.038			0.026			0.028

\* When using a 0.2 L/min sampling rate for 10 min and comparing results to the Transitional Ceiling PEL of 20 ppm.

CV1 (pooled) = 0.031

Withdrawn  
Provided for Historical Reference only

Table 2  
Sampling and Analysis - Hydrogen Sulfide

Test Level*	-----Found-----				Taken ppm	Percent Recovery
	µg	L Air	mg/m3	ppm		
0.5 × PEL	22.33	2.047	10.9	9.30	10.06	92.4
	21.54	2.013	10.7	9.12	10.06	90.7
	11.07	1.967	5.63	4.80	10.06	47.7**
	21.49	1.955	11.0	9.37	10.06	93.1
	21.69	2.007	10.8	9.21	10.06	91.6
	23.75	2.019	11.8	10.03	10.06	99.7
		N	5			
		Mean	9.41			
		Std Dev	0.36			93.5
		CV2	0.038			
1 × PEL	51.10	2.043	25.01	21.38	22.18	96.4
	49.41	2.028	24.36	20.83	22.18	93.9
	29.61	1.956	15.14	12.94	22.18	58.3**
	49.27	1.966	25.06	21.43	22.18	96.6
	51.89	1.993	26.04	22.26	22.18	100.4
	55.36	2.037	27.18	23.23	22.18	104.7
		N	5			
		Mean	21.83			
		Std Dev	0.94			98.4
		CV2	0.043			
2 × PEL	93.00	2.098	44.33	37.64	39.22	96.0
	93.86	2.056	45.65	38.77	39.22	98.8
	91.04	2.013	45.23	38.41	39.22	97.9
	90.07	2.013	44.74	38.00	39.22	96.9
	94.62	2.055	46.04	39.10	39.22	99.7
	99.67	2.083	47.85	40.63	39.22	103.6
		N	6			
		Mean	38.76			
		Std Dev	1.06			98.8
		CV2	0.027			

\* 0.2 L/min sampling rate (10 min sample) and 20 ppm Transitional Ceiling PEL

\*\* Deleted as an outlier (99% confidence level)

CV2 (pooled) = 0.036, CVT (pooled) = 0.038, Average Recovery = 96.9%



Withdrawn  
 Provided for Historical Reference only

Table 3  
 Collection Efficiency - Hydrogen Sulfide (85 to 88% RH)

Sample	ppm H <sub>2</sub> S In Filter	ppm H <sub>2</sub> S In Backup Impinger	Collection Efficiency (%)
1	37.64	ND	100
2	38.77	ND	100
3	38.41	ND	100
4	38.00	ND	100
5	39.10	ND	100
6	40.63	ND	100

ND = No H<sub>2</sub>S detected, detection limit = 0.4 ppm

Table 4  
 Breakthrough - Hydrogen Sulfide (18% RH)

Generated Conc (ppm)	Sampling Time (min)	ppm Found In Backup Impinger		Breakthrough (%)
			CV	
39.78	10	ND	----	0.0
39.78	15	0.98	0.051	2.5
39.78	20	2.22	0.096	5.6

Breakthrough - Hydrogen Sulfide (86% RH)

Generated Conc (ppm)	Sampling Time (min)	ppm Found In Backup Impinger		Breakthrough (%)
			CV	
39.04	10	ND	----	0.0
39.04	15	1.10	0.103	2.8
39.04	20	2.48	0.106	6.3

ND = No H<sub>2</sub>S detected, detection limit = 0.4 ppm

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
 Provided for Historical Reference only

Table 5  
 Storage Stability Test - Hydrogen Sulfide

Storage Day	-----Found-----				Taken ppm	Percent Recovery
	µg	L Air	mg/m3	ppm		
0	51.10	2.043	25.01	21.38	22.18	96.4
	49.41	2.028	24.36	20.83	22.18	93.9
	29.61	1.956	15.14*	12.94*	22.18	58.3*
	49.27	1.966	25.06	21.43	22.18	96.6
	51.89	1.993	26.04	22.26	22.18	100.4
	55.36	2.037	27.18	23.23	22.18	104.7
			N Mean Std Dev CV	5 21.83 0.94 0.043		98.4
5	44.00	1.960	22.45	19.32	20.28	95.3
	49.63	2.114	23.48	20.20	20.28	99.6
	45.92	1.987	23.11	19.88	20.28	98.0
	47.76	2.073	23.04	19.83	20.28	97.8
	47.51	2.070	22.95	19.75	20.28	97.4
	46.88	2.041	22.97	19.76	20.28	97.4
			N Mean Std Dev CV	6 19.79 0.28 0.014		98.4
14	45.46	1.947	23.35	19.98	20.30	98.4
	48.65	2.099	23.18	19.83	20.30	97.7
	45.56	1.971	23.12	19.78	20.30	97.4
	47.55	2.073	22.94	19.63	20.30	96.7
	49.17	2.051	23.97	20.51	20.30	101.0
	48.92	2.026	24.15	20.66	20.30	101.8
			N Mean Std Dev CV	6 20.07 0.42 0.021		98.8
30	49.24	2.028	24.28	20.89	20.28	103.0
	52.14	2.177	23.95	20.60	20.28	101.6
	49.29	2.042	24.14	20.76	20.28	102.4
	49.82	2.144	23.24	19.99	20.28	98.6
	50.70	2.134	23.76	20.44	20.28	100.8
	50.78	2.104	24.14	20.76	20.28	102.4
			N Mean Std Dev CV	6 20.57 0.32 0.016		101.4

\* Deleted as an outlier.

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
 Provided for Historical Reference only

Table 6  
 Quantitative Detection Limit - Hydrogen Sulfide

-----Amount of Hydrogen Sulfide-----						
Sample No.	-----10.2 ng/mL-----		-----20.4 ng/mL-----		-----50.9 ng/mL-----	
	ng/mL found	response (μA)	ng/mL found	response (μA)	ng/mL found	response (μA)
1	8.5	0.0198	19.5	0.0527	51.8	0.1473
2	7.3	0.0159	24.2	0.0669	51.2	0.1455
3	-	ND	22.2	0.0607	54.1	0.1538
4	9.5	0.0228	20.3	0.0553	51.9	0.1475
5	-	ND	17.9	0.0481	53.3	0.1515
ND = None detected						
N	5	5	5	5	5	5
Mean	-	-	20.8	0.0567	52.5	0.1491
Std Dev	-	-	2.4	0.0073	1.2	0.0034
CV	-	-	0.117	0.128	0.023	0.023

-----Amount of Hydrogen Sulfide-----				
Sample No.	-----101.8 ng/mL-----		-----203.5 ng/mL-----	
	ng/mL found	response (μA)	ng/mL found	response (μA)
1	100.6	0.2814	206.0	0.5370
2	100.4	0.2810	207.4	0.5400
3	102.2	0.2857	212.4	0.5510
4	98.9	0.2769	198.8	0.5210
5	100.5	0.2812	194.3	0.5110
N	5	5	5	5
Mean	100.5	0.2812	203.9	0.5320
Std Dev	1.7	0.0031	7.2	0.0159
CV	0.012	0.011	0.035	0.030

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
 Provided for Historical Reference only

Table 7  
 Independent Laboratory Assessment of Method at Low Humidity  
 Hydrogen Sulfide

Sample #	-----Found-----				Taken ppm	Percent Recovery
	µg	L Air	mg/m3	ppm		
5	16.36	2	8.18*	5.78*	10	57.9
6	36.72	2	18.36	12.98	10	129.8
7	34.12	2	17.06	12.06	10	120.6
8	31.86	2	15.93	11.26	10	112.6
9	33.72	2	16.86	11.92	10	119.2
10	31.81	2	15.90	11.24	10	112.4
			N	5		
			Mean	11.89		118.9
			Std Dev	0.71		7.1
			CV	0.060		
11	59.15	2	29.58	20.91	20	104.6
12	58.23	2	29.12	20.58	20	102.9
13	58.23	2	29.12	20.58	20	102.9
14	55.87	2	27.94	19.75	20	98.8
15	56.71	2	28.36	20.05	20	100.2
16	55.18	2	27.59	19.50	20	97.5
			N	6		
			Mean	20.23		101.1
			Std Dev	0.55		2.7
			CV	0.027		

Results for six samples taken at 40 ppm H<sub>2</sub>S and under the same conditions gave recoveries about 50% lower than theoretical; however, generation system problems were discovered and the lower recovery was attributed to a faulty H<sub>2</sub>S permeation tube used to calibrate the GC. The GC was used to monitor and determine the known concentration of the generation.

\*Deleted as an outlier.

Withdrawn  
 Provided for Historical Reference only

These samples were taken at atmospheric conditions of 761 mmHg and 21 °C.

Sample #	-----Found-----				Taken ppm	Percent Recovery
	µg	L Air	mg/m3	ppm		
1	32.2	2	16.1	11.4	12	95.0
2	29.4	2	14.7	10.4	12	86.7
3	37.5	2	18.8	13.3	12	110.8
4	35.7	2	17.8	12.6	12	105.0
5	38.5	2	19.2	13.6	12	113.3
6	36.3	2	18.2	12.9	12	107.5
			N	6		
			Mean	12.4		103.3
			Std Dev	1.2		10.1
			CV	0.098		
7	91.5	2	45.8	32.4	30	108.0
8	93.1	2	46.6	32.9	30	109.7
9	85.5	2	42.8	30.2	30	100.7
10	88.6	2	44.3	31.3	30	104.3
11	97.0	2	48.5	34.3	30	114.3
12	86.3	2	43.2	30.5	30	101.7
			N	3		
			Mean	31.9		106.3
			Std Dev	1.6		5.2
			CV	0.048		

Results for three samples taken at about 38 ppm H<sub>2</sub>S and under the same conditions gave recoveries about 50% lower than theoretical; however, the sampling media had been incorrectly handled prior to sample collection and the lower recovery was attributed to improper handling.

These samples were taken at atmospheric conditions of 761 mmHg and 21 °C.

Appendix

Generation System for Producing Dynamic Test Atmospheres of Hydrogen Sulfide

A generation system was designed such that the H<sub>2</sub>S gas and diluent air were connected to a mixing tee which was then connected to a sampling manifold. The diluent air was conditioned to the RH and temperature required for the particular experiment. Details of the system are listed below.

A. Description of equipment used to monitor and control gases.

1. The flow rate, temperature, and RH of the diluent air were controlled with a model HCS-201 Mass Flow, Temperature, and Humidity Control System (Miller-Nelson Research Inc., Monterey, CA). The humidity and temperature of the air were measured with a model 400D % Relative Humidity/Temperature Monitor (General Eastern, Watertown, MA).
2. The flow rate of the H<sub>2</sub>S gas was controlled with mass flow controllers (Tylan Corp., Carson, CA).

Note: OSHA no longer uses or supports this method (2021)

Withdrawn  
Provided for Historical Reference only

3. The flow rate of the H<sub>2</sub>S gas and the air were measured just prior to and just after sampling. The flow rate of the H<sub>2</sub>S gas was measured with a model 823-1 Electronic Bubble Flowmeter (Mast Development Co., Davenport, IA). The flow rate of the humid air was measured with a dry test gas meter

B. Description of sampling manifold and connective tubing:

1. A 12-port sampling manifold was constructed from Teflon tubing and fittings. There were six sampling ports located on opposite sides of the sampling manifold.
2. The sampling manifold was directly connected to a Teflon mixing tee by means of Teflon tubing and Teflon L fittings.
3. The sampling devices were attached to the manifold with Teflon tubing and to the sampling pumps with Tygon tubing.
4. All other connective tubing and fittings were made of stainless steel or Teflon to avoid any contamination problems.