



## Monocrotophos

---

Method number: PV2045

Target Concentration: 0.25 mg/m<sup>3</sup> OSHA permissible exposure level (PEL).

Procedure: Samples are collected by drawing known volumes of air through OSHA versatile sampler (OVS-2) tubes, each containing a glass fiber filter and two sections of XAD-2 adsorbent. Samples are desorbed with toluene and analyzed by gas chromatography (GC) using a flame photometric detector (FPD).

Recommended air volume and sampling rate: 480 minutes at 1.0 L/min (480 L)

Detection limit of the overall procedure: 7.4 µg/m<sup>3</sup> (based on the recommended air volume and the analytical detection limit):

Status of method: Partially Validated method. This method has been partially evaluated and is presented for information and trial use only.

May 1990 (final)

Duane Lee

Carcinogen and Pesticide Branch  
OSHA Salt Lake Technical Center  
Salt Lake City UT 84115-1802

## 1 General Discussion

### 1.1 Background

#### 1.1.1 History of procedure

This evaluation was undertaken because OSHA recently adopted the TLV's as PEL's. The OVS-2 sampler was tested as an effective sampling device for monocrotophos. This method follows the procedures developed for other organophosphorus pesticides. (Ref. 5.1)

#### 1.1.2 Toxic effects (This section is for information only and should not be taken as the basis of OSHA policy.)

Monocrotophos is an organophosphorus pesticide. It has an oral LD<sub>50</sub> of 17 mg/kg in male rats and 20 mg/kg in female rats. (Ref. 5.4) It is a highly toxic, direct acting, cholinesterase inhibitor, which appears to be capable of penetration through the intact skin; but is excreted rapidly and does not appear to accumulate within the body (Ref. 5.5).


#### 1.1.3 Potential workplace exposure

There was no information available on the number of workers exposed to monocrotophos which is used as a systemic insecticide. There was an estimated  $2.4 \times 10^9$  grams of monocrotophos produced in the U.S. in 1974. (Ref. 5.6)

#### 1.1.4 Physical properties (Ref. 5.2 to 5.5)

CAS number: 6923-22-4  
IMIS number: 2690  
Molecular weight: 233.16  
Molecular formula: C<sub>7</sub>H<sub>14</sub>NO<sub>5</sub>P  
Melting point: 54–55 °C  
Boiling point: 125 °C  
Solubility: soluble in ethanol, acetone, and water. Practically insoluble in diesel oils and kerosene  
Chemical name: dimethyl 2-methylcarbamoyl-1-methyl-vinyl phosphate  
Synonyms: Apadrin, Azodrin, Biloborn, C1414, Crisodrin, Hazodrin, SO9129, Monocron, Nuvacron, Plantdrin, Pilladrin  
Description: reddish brown solid with a mild ester odor

Structure:



### 1.2 Limit defining parameters

The detection limit of the analytical procedure, including a 28:1 split ratio, is 0.063 ng per injection. This is the amount of analyte which will give a peak whose height is approximately five times the baseline noise. (Figure 1.)

## 2 Sampling Procedure

### 2.1 Apparatus

- 2.1.1 A personal sampling pump that can be calibrated to within  $\pm 5\%$  of the recommended flow rate with the sampling device in line.
- 2.1.2 OVS-2 tubes, which are specially made 13-mm o.d. glass tubes that are tapered to 6-mm o.d. They are packed with two sections of cleaned XAD-2 adsorbent and a 13-mm diameter glass fiber filter. The sampling section and backup section contain 270 and 140 mg adsorbent respectively. The backup section is retained by two foam plugs, and the sampling section is between a foam plug and the glass fiber filter. The glass fiber filter is held next to the sampling section by a polytetrafluoroethylene (PTFE) retainer. (Figure 2.)

### 2.2 Reagents

No sampling reagents are required.

### 2.3 Sampling technique

- 2.3.1 Immediately before sampling, remove the plastic caps from the OVS-2 tube.
- 2.3.2 Attach the small end of the tube to the sampling pump with flexible tubing.
- 2.3.3 Attach the tube vertically in the employee's breathing zone in such a manner that it does not impede work performance.
- 2.3.4 After sampling for the appropriate time, remove the tube and seal it with the plastic caps.
- 2.3.5 Wrap each sample end-to-end with a Form OSHA-21 seal.
- 2.3.6 Record the air volume for each sample, and list any possible interference.
- 2.3.7 Submit at least one blank for each set of samples. Handle the blank in the same manner as the samples, except no air is drawn through it.
- 2.3.8 Submit bulk samples for analysis in a separate container. Do not ship with air samples.

### 2.4 Desorption efficiency (glass fiber filter and XAD-2 adsorbent)

Six 4-mL vials each containing a 13-mm glass fiber filter and 270-mg of XAD-2 adsorbent were each liquid spiked on the glass fiber filter with 8.2  $\mu\text{L}$  of a 14.7 mg/mL standard of monocrotophos and allowed to dry overnight in a drawer at ambient temperature. These samples were each desorbed with 2.0 mL of toluene containing triphenyl phosphate (TPP) as the internal standard, shaken for 30 min, and analyzed as in Section 3. The results are listed in Table 2.4.

Table 2.4  
Desorption Efficiency

sample #	µg spiked	µg found	% recovered
Ex1	120.54	112.66	93.5
Ex2	120.54	117.83	97.7
Ex3	120.54	107.09	88.8
Ex4	120.54	118.76	98.5
Ex5	120.54	117.59	97.5
Ex6	120.54	116.49	96.6

average = 95.4%

## 2.5 Retention efficiency

Eighteen OVS-2 tubes were each liquid spiked with 8.2 µL of a 14.7 mg/mL monocrotophos standard on the glass fiber filter. They were allowed to dry overnight and then 480 L of humid air (~77% relative humidity) was drawn through each tube at 1 L/min. Six of the tubes were each desorbed with 2.0 mL of toluene containing TPP, shaken for 30 min and then analyzed as in Section 3. The results are listed in Table 2.5. The remaining samples were stored, 6 in a drawer at ambient temperature and 6 in a freezer.

Table 2.5  
Retention Efficiency

sample #	µg spiked	µg found	% recovered
R1	120.54	99.75	82.8
R2	120.54	110.41	91.6
R3	120.54	96.83	80.3
R4	120.54	96.06	79.7
R5	120.54	96.67	80.2
R6	120.54	106.7	88.5

average = 83.8%

## 2.6 Sample storage

After 3 days of storage, 6 tubes, 3 from the ambient storage group and 3 from the freezer storage group, were each desorbed with 2.0 mL of toluene containing TPP, shaken for 30 min and then analyzed as in Section 3. The remaining tubes were desorbed and analyzed after 9 or 10 days of storage. The results are given in Tables 2.6.1 and 2.6.2.

Table 2.6.1  
Ambient Storage

sample #	µg spiked	µg found	% recovered
3	120.54	99.8	82.8
3	120.54	100.59	83.4
4	120.54	93.2	77.3
9	120.54	92.42	76.7
10	120.54	108.03	89.6
10	120.54	99.43	82.5

Average of ~3 days = 81.2%

Average of ~10 days = 82.9%

Table 2.6.2  
Freezer Storage

sample #	µg spiked	µg found	% recovered
3	120.54	101.16	83.9
3	120.54	104.47	86.7
3	120.54	104.09	86.4
9	120.54	102.14	84.7
9	120.54	104.0	86.3
9	120.54	105.14	87.2

Average of 3 days = 85.7%

Average of 9 days = 86.1%

## 2.7 Recommended air volume and sampling rate

2.7.1 The recommended air volume is 480 L.

2.7.2 The recommended flow rate is 1.0 L/min.

## 2.8 Interferences (sampling)

It is not known if any compounds will interfere with the collection of monocrotophos. Any suspected interferences should be reported to the laboratory.

## 2.9 Safety precautions (sampling)

2.9.1 Attach the sampling equipment in such a manner that it will not interfere with work performance or employee safety.

2.9.2 Follow all safety practices that apply to the work area being sampled.

## 3 Analytical Procedure

### 3.1 Apparatus

3.1.1 A balance capable of weighing to the nearest tenth of a milligram. A Mettler HL52 balance was used in this evaluation.

3.1.2 A mechanical shaker.

3.1.3 A GC equipped with an FPD. A Hewlett Packard (HP) 5890, equipped with an FPD and an autosampler, was used in this evaluation.

3.1.4 A GC column capable of separating monocrotophos from any interference. A 30-m × 0.32-mm i.d. (1.0 µm d<sub>r</sub> DB-5) capillary column was used in this evaluation.

3.1.5 An electronic integrator or some other suitable means for measuring detector response. A Waters 860 Laboratory Data System was used in this evaluation.

3.1.6 Volumetric flasks and pipettes.

3.1.7 Vials, 2-mL and 4-mL.

### 3.2 Reagents

3.2.1 Toluene, reagent grade.

3.2.2 Monocrotophos, reagent grade. A standard obtained from EPA (EPA #0360, 98% purity) was used in this evaluation.

3.2.3 Triphenyl phosphate (TPP), reagent grade. A 40 µg/mL solution in toluene was used as an internal standard (Istd).

### 3.3 Standard preparation

Prepare monocrotophos stock standards by weighing 10 to 15 mg of monocrotophos. Transfer the monocrotophos to separate 10-mL volumetric flasks, and add toluene to the mark. Make working range standards of 1.7 to 147 µg/mL by pipette dilutions of the stock standards with toluene containing TPP. Store stock and dilute standards in a freezer.

### 3.4 Sample preparation

3.4.1 Transfer the 13-mm glass fiber filter and the 270-mg sampling section of the tube to a 4-mL vial. Place the first foam plug and the 140-mg section in a separate 4-mL vial. A small glass funnel can be used to facilitate the transfer of the adsorbent. Discard the rear foam plug. Do not discard the glass sampling tube; it can be reused.

3.4.2 Add 2.0 mL of toluene containing TPP to each vial and seal with a PTFE-lined cap.

3.4.3 Shake the vials for 30 minutes on a mechanical shaker.

3.4.4 Transfer, if necessary, the samples to 2-mL vials for use on an HP autosampler.

### 3.5 Analysis

#### 3.5.1 Instrument conditions

**Column:** 30 m × 0.32 mm i.d., (1.0 µm d<sub>f</sub> DB-5)

**Temperatures:**

Injector temperature: 250 °C

Column temperature: 220 °C

Detector temperature: 225 °C

**Gas flows:** mL/min

Column (Hydrogen): 1

FPD make up (nitrogen): 42

**Injection volume:** 1 µL

**Split ratio:** 28:1

**Retention time:** 4.66 min

#### 3.5.2 Chromatogram (Figure 3.)

### 3.6 Interferences (analytical)

3.6.1 Any collected compound having a similar retention time to that of the analyte or the internal standard is a potential interference.

3.6.2 GC conditions may generally be varied to circumvent interferences.

3.6.3 Retention time on a single column is not proof of chemical identity. Analysis by an alternate GC column, high performance liquid chromatography (HPLC), and confirmation by mass spectrometry are additional means of identification.

### 3.7 Calculations

3.7.1 Construct a calibration curve (Figure 4.) by plotting detector response versus concentration ( $\mu\text{g}/\text{mL}$ ) of monocrotophos.

3.7.2 Determine the  $\mu\text{g}/\text{mL}$  of monocrotophos in both sections of each sample and blank from the calibration curve.

3.7.3 Blank correct each section by subtracting the  $\mu\text{g}/\text{mL}$  found in each blank section from the  $\mu\text{g}/\text{mL}$  found in each corresponding sample section and then add the values together.

3.7.4 Determine the air concentration by using the following formula.

$$\text{mg} / \text{m}^3 = \frac{(\mu\text{g} / \text{mL, blank corrected})(\text{desorption volume, mL})}{(\text{air volume, L})(\text{desorption efficiency, decimal})}$$

### 3.8 Safety precautions (analytical)

3.8.1 Avoid skin contact and air exposure to monocrotophos.

3.8.2 Avoid skin contact with all solvents.

3.8.3 Wear safety glasses at all times.

## 4 Recommendation for Further Study

This method should be fully validated.




Figure 1.  
Detection Limit Chromatogram for Monocrotophos, Crufomate, and TPP (ISTD)




Figure 2.  
OVS-2 Sampling Tube




Figure 3.  
Chromatogram of Monocrotophos with Crufomate and TPP (Istd)




Figure 4.  
Calibration Curve for Monocrotophos



## 5 References

- 5.1 Burreight, D., Method #62, Chlorpyrifos, DDVP, Diazinon, Malathion, and Parathion, OSHA Analytical Laboratory, unpublished, 1986.
- 5.2 Registry of Toxic Effects of Chemical Substances 1985-86 Edition; DHHS (NIOSH) Publication No. 87-114, U.S. Department of Health and Human Services: Cincinnati, OH, 1987; p 3392.
- 5.3 Farm Chemicals Handbook; Berg, Gordon L. Ed.; Meister: Willoughby, Ohio, 1989; p C201.
- 5.4 Merck Index, 10th ed.; Windholz, Martha ED.; Merck: Rahway, N.J., 1983; p 894.
- 5.5 Documentation of Threshold Limit Values and Biological Exposure Indices; American Conference of Governmental Industrial Hygienists Inc., Fifth Edition, 1986, p 416.
- 5.6 HSDB (Hazardous Substance Data Base).